
colloidal gold-Rabbit anti Goat IgG(H+L)

Product Information

Catalog No.:	FNSA-0102
Size:	100µg/200µl
Form:	liquid
Host:	Rabbit
Purification:	The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Fluorophore:	colloidal gold(40nm)
Storage:	2mM boric acid, 1%BSA, 1%trehalose, 0.1%NaN ₃ , pH9.0, 4°C for 24 months
Dilution buffer:	2mM boric acid, 1%BSA, 0.1%NaN ₃ , pH9.0 ProClin300 should be avoided.

Note: This product is for research use only, not for diagnostic or therapeutic use.

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone, Wuhan, Hubei, China(430206)

Tel:(0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com