

---

## HRP conjugated Mouse Anti-Rabbit IgG Light Chain Specific

### Product Information

Catalog No.:	FNSA-0110
Size:	100µl/500µl
Form:	liquid
Purification:	The antibody was purified by chromatographic procedures.
Conjugation:	HRP purified from horseradish root
Stabilizer:	3mg/ml BSA
Storage:	0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6. -20°C for 24 months.

### Application

Tested Application: ELISA, WB

Suggested dilution: 1:4000-1:10000 for western blotting with ECL substrates

1:2000-1:50000 for ELISA and Western blotting with chromogenic substrates

**Note: This product is for research use only, not for diagnostic or therapeutic use.**