

HRP conjugated Mouse Anti-Rabbit IgG Heavy Chain Specific

Product Information

Catalog No.:	FNSA-0111
Size:	100μl/500μl
Form:	liquid
Purification:	The antibody was affinity purified.
Conjugation:	HRP purified from horseradish root
Stabilizer:	0.5% BSA
Storage:	0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, 0.1% proclin300, pH7.3. -20°C for 24 months.

Application

Tested Application: ELISA, WB

Suggested dilution: 1:5000-1:20000 for western blotting with ECL substrates

1:5000-1:20000 for ELISA and Western blotting with chromogenic substrates

Note: This product is for research use only, not for diagnostic or therapeutic use.

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com