

---

## **Biotin conjugated 2019 nCOV N protein antibody (1G5)**

### **Product Information**

Catalog No.:	FNSA-0115
Size:	100µL/500µL
Host:	Mouse
Isotype:	IgG1
Form:	liquid
Conjugation:	Biotin
Stabilizer:	5mg/ml BSA ( IgG-Free, Protease-Free)
Storage:	0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6. -20℃ for 24 months , avoid repeated freezing and thawing.

### **Application**

Tested Application:ELISA, WB, IHC

**Note: This product is for research use only, not for diagnostic or therapeutic use.**