

## Goat Anti-Human Whole Serum

### Product Information

Catalog No.:	FNSA-0117
Size:	1mg
Form:	liquid
Source:	Goat
Clone:	polyclonal
Stabilizer:	3mg/ml BSA ( IgG-Free, Protease-Free)
Storage:	0.01M phosphate buffered saline, pH7.4, containing 0.1% proclin300 as preservative. -20°C for 24 months, avoid repeated freezing and thawing.

### Application

Tested Application: ELISA, WB, IHC, IP

**Note: This product is for research use only, not for diagnostic or therapeutic use.**