
Normal Mouse IgG Control

Product Information

Catalog No.:	FNSA-0130
Size:	100μl/ 500μl
Form:	liquid
Source:	mouse serum
Purification:	Protein A +G purified
Storage:	0.01M Sodium phosphate, 0.25M NaCl, 0.1% proclin300, 50% glycerol, pH7.4. -20℃ for 24 months , avoid repeated freezing and thawing.

Application

Negative control for use in conjunction with other antibodies. Use at the same concentration as the detection antibody.

Note: This product is for research use only, not for diagnostic or therapeutic use.

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone, Wuhan, Hubei, China (430206)

Tel: (0086) 027-87384275

Fax: (0086) 027-87800889

www.fn-test.com