

## Chicken IgY- Isotype Control

### Product Information

Catalog No.:	FNSA-0131
Size:	1mg
Form:	Liquid
Clone:	polyclonal
Isotype:	IgY
Purity:	≥95% as determined by SDS-PAGE
Source:	Normal Chicken IgY is obtained from chickens that have not been Immunized.
Storage:	0.01M Sodium phosphate, 0.25M NaCl, 0.1% proclin300, pH 7.4. 12 months from date of receipt at -4℃.

### Application

Negative control for use in conjunction with other antibodies. Use at the same concentration as the detection antibody.

**Note: This product is for research use only, not for diagnostic or therapeutic use.**

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)