
Alexa Fluor 555-Goat Anti-Mouse IgG(H+L)

Product Information

Catalog No.:	FNSA-0156
Size:	100µl/500µl
Form:	liquid
Purification:	The antibody was purified by immunoaffinity chromatography
Conjugation:	Alexa Fluor 555
Fluorophore:	Alexa Fluor 555 carboxylic acid
Host:	Goat
Storage:	0.01M Sodium phosphate, 0.25M NaCl, 0.02% sodium azide, 50% glycerol, pH7.6. -20°C for 24 months, avoid repeated freezing and thawing. Store in dark.

Application

Tested Application:	IF, FC
Recommended dilution:	1:100-1:1000

Note: This product is for research use only, not for diagnostic or therapeutic use.