
HRP-Goat Anti-Human IgG (Fab)

Product Information

| | |
|--------------|--|
| Catalog No.: | FNSA-0160 |
| Size: | 100µl/500µl |
| Form: | liquid |
| Conjugation: | HRP purified from horseradish root |
| Stabilizer: | 3mg/ml BSA (IgG-Free, Protease-Free) |
| Storage: | 0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6. -20°C for 24 months , avoid repeated freezing and thawing. |

Application

Tested Application:ELISA, WB

Suggested dilution: 1:2000-1:10000 for western blotting with ECL substrates

1:1000-1:20000 for ELISA and Western blotting with chromogenic substrates

Note: This product is for research use only, not for diagnostic or therapeutic use.