

---

## HRP conjugated Goat Anti-Human IgG (Fab)

### Product Information

|              |  |
|--------------|--|
| Catalog No.: | FNSA-0160  |
| Size:        | 100µl/500µl  |
| Form:        | liquid   |
| Conjugation: | HRP purified from horseradish root   |
| Stabilizer:  | 3mg/ml BSA ( IgG-Free, Protease-Free)  |
| Storage:     | 0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6.<br>-20°C for 24 months , avoid repeated freezing and thawing. |

### Application

Tested Application:ELISA, WB

Suggested dilution: 1:2000-1:10000 for western blotting with ECL substrates

1:1000-1:20000 for ELISA and Western blotting with chromogenic substrates

**Note: This product is for research use only, not for diagnostic or therapeutic use.**