

---

## Goat anti- human IgG(Fab)

### Product Information

Catalog No.:	FNSA-0164
Size:	1mg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Goat
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide , -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen: human IgG(Fab)

### Application

Reactivity: Human IgG(Fab)