
Biotin conjugated anti-Mouse IP10

Product Information

Catalog No.:	FNSA-0169
Size:	100µL/500µL
Host:	Rabbit
Clonality:	polyclonal
Isotype:	IgG
Form:	liquid
Conjugation:	Biotin
Stabilizer:	5mg/ml BSA (IgG-Free, Protease-Free)
Storage:	0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6. -20°C for 24 months , avoid repeated freezing and thawing.

Application

Tested Application:ELISA

Note: This product is for research use only, not for diagnostic or therapeutic use.