

---

## Alexa Fluor 647-Goat Anti-Rabbit IgG(H+L)

### Product Information

Catalog No.:	FNSA-0177
Size:	100µl/500µl
Form:	liquid
Purification:	The antibody was purified by immunoaffinity chromatography
Conjugation:	Alexa Fluor 647
Fluorophore:	Alexa Fluor 647, Amax=654 nm, Emax=674 nm
Host:	Goat
Storage:	PBS with 50% glycerol, 5 mg/mL BSA, 0.1% Proclin-300, pH 7.4. -20°C for 24 months, avoid repeated freezing and thawing. Store in dark.

### Application

Tested Application:	IF, FC
Recommended dilution:	1:100-1:1000

**Note: This product is for research use only, not for diagnostic or therapeutic use.**