
Poly-HRP-Mouse Anti-Human IgA

Product Information

Catalog No.:	FNSA-0178
Size:	100µl/500µl
Form:	liquid
Purification:	The antibody was purified by immunoaffinity chromatography
Conjugation:	Poly-HRP
Host:	Mouse
Stabilizer:	3mg/ml BSA (IgG-Free, Protease-Free)
Storage:	0.01M Sodium phosphate, 0.25M NaCl, 0.1% proclin300, 50% glycerol, pH7.4. -20 °C for 24 months, avoid repeated freezing and thawing.
Application:	ELISA
Recommended dilution:	1:5000-1:20000

Note: This product is for research use only, not for diagnostic or therapeutic use.