

## A1BG antibody

### Product Information


|               |  |
|---------------|--|
| Catalog No.:  | FNab00012  |
| Size:         | 100µg  |
| Form:         | liquid   |
| Purification: | Immunogen affinity purified  |
| Purity:       | ≥95% as determined by SDS-PAGE   |
| Host:         | Rabbit   |
| Clonality:    | polyclonal   |
| Clone ID:     | None   |
| IsoType:      | IgG  |
| Storage:      | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.) |

### Immunogen information


|              |   |
|--------------|---|
| Immunogen:   | alpha-1-B glycoprotein                            |
| Synonyms:    | Alpha-1B-glycoprotein Alpha-1-B glycoprotein A1BG |
| Observed MW: | 65 kDa  |
| Uniprot ID : | P04217  |

### Application

|                       |   |
|-----------------------|---|
| Reactivity:           | Human, Mouse  |
| Tested Application:   | ELISA, WB, IHC, IF                                  |
| Recommended dilution: | WB: 1:2000-1:20000; IF: 1:10-1:100; IHC: 1:10-1:200 |
| Image:                |   |



Immunohistochemistry of paraffin-embedded human hepatocirrhosis using FNab00012(A1BG antibody) at dilution of 1:50



human plasma tissue were subjected to SDS PAGE followed by western blot with FNab00012(A1BG Antibody) at dilution of 1:5000