

## A4GNT antibody

### Product Information

Catalog No.:	FNab00015
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Necessary for the synthesis of type III mucin. Catalyzes the transfer of N-acetylglucosamine(GlcNAc) to core 2 branched O-glycans.

### Immunogen information


Immunogen:	alpha-1,4-N-acetylglucosaminyltransferase
Synonyms:	Alpha-1,4-N-acetylglucosaminyltransferase (Alpha4GnT) A4GNT
Observed MW:	39 kDa
Uniprot ID :	Q9UNA3

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:5000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human stomach slide using FNab00015(A4GNT Antibody) at dilution of 1:50



K-562 cells were subjected to SDS PAGE followed by western blot with FNab00015(A4GNT antibody) at dilution of 1:1000