

AACS antibody

Product Information

Catalog No.:	FNab00016
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Activates acetoacetate to acetoacetyl-CoA. May be involved in utilizing ketone body for the fatty acid-synthesis during adipose tissue development(By similarity).

Immunogen information

Immunogen:	acetoacetyl-CoA synthetase
Synonyms:	Acetoacetyl-CoA synthetase Acyl-CoA synthetase family member 1 Protein sur-5 homolog AACS ACSF1
Observed MW:	76 kDa
Uniprot ID :	Q86V21

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:200-1:1000; IF: 1:20-1:200
Image:	



mouse heart tissue were subjected to SDS PAGE followed by western blot with FNab00016(AACS Antibody) at dilution of 1:300