

ABCA8 antibody

Product Information

Catalog No.:	FNab00029
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


ATP-dependent lipophilic drug transporter.

Immunogen information

Immunogen:	ATP-binding cassette, sub-family A(ABC1), member 8
Synonyms:	ABC-type organic anion transporter ABCA8 ATP-binding cassette sub-family A member 8 ABCA8 KIAA0822
Observed MW:	179 kDa
Uniprot ID :	O94911

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:200-1:1000
Image:	



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab00029(ABCA8 Antibody) at dilution of 1:300