

## ABCB6 antibody

### Product Information

Catalog No.:	FNab00032
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Binds heme and porphyrins and functions in their ATP-dependent uptake into the mitochondria. Plays a crucial role in heme synthesis.

### Immunogen information


Immunogen:	ATP-binding cassette, sub-family B(MDR/TAP), member 6
Synonyms:	ATP-binding cassette sub-family B member 6 ABC-type heme transporter ABCB6 Mitochondrial ABC transporter 3 (Mt-ABC transporter 3) P-glycoprotein-related protein Ubiquitously-expressed mammalian ABC half transporter ABCB6 MTABC3 PRP UMAT
Observed MW:	94 kDa
Uniprot ID :	Q9NP58

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC,IF
Recommended dilution:	WB: 1:500-1:5000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer using FNaB00032(ABCB6 antibody) at dilution of 1:50



Raji cells were subjected to SDS PAGE followed by western blot with FNaB00032(ABCB6 antibody) at dilution of 1:500