

ABCD1 antibody

Product Information

Catalog No.:	FNab00036
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Probable transporter. The nucleotide-binding fold acts as an ATP-binding subunit with ATPase activity.

Immunogen information

Immunogen:	ATP-binding cassette, sub-family D(ALD), member 1
Synonyms:	ATP-binding cassette sub-family D member 1 Adrenoleukodystrophy protein (ALDP) ABCD1 ALD
Observed MW:	80 kDa
Uniprot ID :	P33897

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:5000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab00036(ABCD1 antibody) at dilution of 1:1000