

ABHD1 antibody

Product Information


Catalog No.:	FNab00044
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	abhydrolase domain containing 1
Synonyms:	Protein ABHD1 Alpha/beta hydrolase domain-containing protein 1 (Abhydrolase domain-containing protein 1) Lung alpha/beta hydrolase 1 ABHD1 LABH1
Observed MW:	45 kDa
Uniprot ID :	Q96SE0

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:200-1:2000
Image:	



mouse heart tissue were subjected to SDS PAGE followed by western blot with FNab00044(ABHD1 Antibody) at dilution of 1:600