

## ABHD14B antibody

### Product Information

Catalog No.:	FNab00047
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Has hydrolase activity towards p-nitrophenyl butyrate(in vitro). May activate transcription.

### Immunogen information


Immunogen:	abhydrolase domain containing 14B
Synonyms:	Putative protein-lysine deacylase ABHD14B Alpha/beta hydrolase domain-containing protein 14B (Abhydrolase domain-containing protein 14B) CCG1-interacting factor B ABHD14B CIB
Observed MW:	22 kDa
Uniprot ID :	Q96IU4

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:5000; IHC: 1:20-1:200; IF: 1:10-1:100
Image:	



Immunohistochemistry of paraffin-embedded human prostate hyperplasia using FNab00047(ABHD14B antibody) at dilution of 1:50



L02 cells were subjected to SDS PAGE followed by western blot with FNab00047(ABHD14B antibody) at dilution of 1:400