

ABLIM1 antibody

Product Information

Catalog No.:	FNab00054
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May act as scaffold protein(By similarity). May play a role in the development of the retina. Has been suggested to play a role in axon guidance.

Immunogen information


Immunogen:	actin binding LIM protein 1
Synonyms:	Actin-binding LIM protein 1 (abLIM-1) Actin-binding LIM protein family member 1 Actin-binding double zinc finger protein LIMAB1 Limatin ABLIM1 ABLIM KIAA0059 LIMAB1
Observed MW:	88 kDa
Uniprot ID :	O14639

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:5000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human liver using FNab00054(ABLIM1 antibody) at dilution of 1:100



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab00054(ABLIM1 antibody) at dilution of 1:1000