

ACAA2 antibody

Product Information

Catalog No.:	FNab00059
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Abolishes BNIP3-mediated apoptosis and mitochondrial damage.

Immunogen information


Immunogen:	acetyl-Coenzyme A acyltransferase 2
Synonyms:	3-ketoacyl-CoA thiolase, mitochondrial Acetyl-CoA acetyltransferase Acetyl-CoA acyltransferase Acyl-CoA hydrolase, mitochondrial Beta-ketothiolase Mitochondrial 3-oxoacyl-CoA thiolase T1 ACAA2
Observed MW:	42 kDa
Uniprot ID :	P42765

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lung cancer using FNA00059(ACAA2 antibody) at dilution of 1:100



HeLa cells were subjected to SDS PAGE followed by western blot with FNA00059(ACAA2 antibody) at dilution of 1:500