

ACAP2 antibody

Product Information

Catalog No.:	FNab00070
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


GTPase-activating protein(GAP) for ADP ribosylation factor 6(ARF6).

Immunogen information

Immunogen:	ArfGAP with coiled-coil, ankyrin repeat and PH domains 2
Synonyms:	Arf-GAP with coiled-coil, ANK repeat and PH domain-containing protein 2 Centaurin-beta-2 (Cnt-b2) ACAP2 CENTB2 KIAA0041
Observed MW:	88 kDa
Uniprot ID :	Q15057

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:200-1:2000; IP: 1:200-1:1000
Image:	



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with FNab00070(ACAP2 antibody) at dilution of 1:300