

ACBD4 antibody

Product Information

Catalog No.:	FNab00075
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Binds medium-and long-chain acyl-CoA esters and may function as an intracellular carrier of acyl-CoA esters.

Immunogen information

Immunogen:	acyl-Coenzyme A binding domain containing 4
Synonyms:	Acyl-CoA-binding domain-containing protein 4 ACBD4
Observed MW:	38 kDa
Uniprot ID :	Q8NC06

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:5000; IF: 1:20-1:200
Image:	



human kidney tissue were subjected to SDS PAGE followed by western blot with FNab00075(ACBD4 Antibody) at dilution of 1:1000