

ACBD5 antibody

Product Information

Catalog No.:	FNab00076
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Acyl-CoA binding protein which acts as the peroxisome receptor for pexophagy but is dispensable for aggrephagy and nonselective autophagy. Binds medium-and long-chain acyl-CoA esters.

Immunogen information

Immunogen:	acyl-Coenzyme A binding domain containing 5
Synonyms:	Acyl-CoA-binding domain-containing protein 5 ACBD5 KIAA1996
Observed MW:	65 kDa
Uniprot ID :	Q5T8D3

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:200-1:2000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab00076(ACBD5 antibody) at dilution of 1:500