

CNNM1 antibody

Product Information

Catalog No.:	FNab00079
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


ACDP1 also known as CNNM1 is restricted to brain and testis. It is probable metal transporter.

Immunogen information

Immunogen:	cyclin M1
Synonyms:	Metal transporter CNNM1 Ancient conserved domain-containing protein 1 Cyclin-M1 CNNM1 ACDP1
Observed MW:	100 kDa
Uniprot ID :	Q9NRU3

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:5000; IP: 1:500-1:5000
Image:	



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab00079(CNNM1 Antibody) at dilution of 1:1000