

ACOT11 antibody

Product Information

Catalog No.:	FNab00089
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Has acyl-CoA thioesterase activity towards medium(C12) and long-chain(C18) fatty acyl-CoA substrates.

Immunogen information

Immunogen:	acyl-CoA thioesterase 11
Synonyms:	Acyl-coenzyme A thioesterase 11 (Acyl-CoA thioesterase 11) Acyl-CoA thioester hydrolase 11 Adipose-associated thioesterase Brown fat-inducible thioesterase (BFIT) Palmitoyl-coenzyme A thioesterase ACOT11 BFIT KIAA0707 THEA
Observed MW:	68.4 kDa
Uniprot ID :	Q8WXI4

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human prostate cancer using FNab00089(ACOT11 antibody) at dilution of 1:100