

ACOX2 antibody

Product Information

Catalog No.:	FNab00096
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Oxidizes the CoA esters of the bile acid intermediates di-and tri-hydroxycholestanoic acids.

Immunogen information

Immunogen:	acyl-Coenzyme A oxidase 2, branched chain
Synonyms:	Peroxisomal acyl-coenzyme A oxidase 2 3-alpha,7-alpha,12-alpha-trihydroxy-5-beta-cholestanoyl-CoA 24-hydroxylase 3-alpha,7-alpha,12-alpha-trihydroxy-5-beta-cholestanoyl-CoA oxidase Trihydroxycoprostanoyl-CoA oxidase (THCA-CoA oxidase, THCCox) ACOX2
Observed MW:	77 kDa
Uniprot ID :	Q99424

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, IF
Recommended dilution:	IF: 1:10-1:100; IHC: 1:20-1:200
Image:	



Immunohistochemical of paraffin-embedded human liver using FNab00096(ACOX2 antibody) at dilution of 1:50