

ACP2 antibody

Product Information


Catalog No.:	FNab00099
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	acid phosphatase 2, lysosomal
Synonyms:	Lysosomal acid phosphatase (LAP) ACP2
Observed MW:	50 kDa, 60 kDa
Uniprot ID :	P11117

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:200-1:2000
Image:	



mouse kidney tissue were subjected to SDS PAGE followed by western blot with FNab00099(ACP2 antibody) at dilution of 1:500