

ACP4 antibody

Product Information

| | |
|---------------|----------------------------------------------------------------------------------------------------------------|
| Catalog No.: | FNab00103 |
| Size: | 100µg |
| Form: | liquid |
| Purification: | Immunogen affinity purified |
| Purity: | ≥95% as determined by SDS-PAGE |
| Host: | Rabbit |
| Clonality: | polyclonal |
| Clone ID: | None |
| IsoType: | IgG |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.) |

Background


Dephosphorylates receptor tyrosine-protein kinase erbB-4 and inhibits the ligand-induced proteolytic cleavage.

Immunogen information

| | |
|--------------|----------------------------------------------------------|
| Immunogen: | acid phosphatase, testicular |
| Synonyms: | Testicular acid phosphatase Acid phosphatase 4 ACP4 ACPT |
| Observed MW: | 35 kDa |
| Uniprot ID : | Q9BZG2 |

Application

| | |
|-----------------------|------------------|
| Reactivity: | Human, Mouse |
| Tested Application: | ELISA, WB |
| Recommended dilution: | WB: 1:500-1:5000 |
| Image: | |



mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab00103(ACPT Antibody) at dilution of 1:1000