
ACSL4/FACL4 antibody

Product Information


Catalog No.:	FNab00108
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	acyl-CoA synthetase long-chain family member 4
Synonyms:	ACS4, ACSL4, FACL4, LACS 4, LACS4, MRX63, MRX68
Observed MW:	79 kDa
Uniprot ID :	O60488

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human placenta slide using FNab00108(ACSL4 Antibody) at dilution of 1:50



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab00108(ACSL4 Antibody) at dilution of 1:1000