

## ACSM5 antibody

### Product Information

Catalog No.:	FNab00112
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Has medium-chain fatty acid:CoA ligase activity with broad substrate specificity(in vitro). Acts on acids from C(4) to C(11) and on the corresponding 3-hydroxy-and 2,3-or 3,4-unsaturated acids(in vitro)(By similarity).

### Immunogen information


Immunogen:	acyl-CoA synthetase medium-chain family member 5
Synonyms:	Acyl-coenzyme A synthetase ACSM5, mitochondrial Acyl-CoA synthetase medium-chain family member 5 ACSM5 MACS3
Observed MW:	65 kDa
Uniprot ID :	Q6NUN0

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using FNab00112(ACSM5 Antibody) at dilution of 1:50



mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab00112(ACSM5 Antibody) at dilution of 1:300