

## ACTL7A antibody

### Product Information

Catalog No.:	FNab00119
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	actin-like 7A
Synonyms:	Actin-like protein 7A Actin-like-7-alpha ACTL7A
Observed MW:	49 kDa
Uniprot ID :	Q9Y615

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:200-1:2000; IF: 1:10-1:100
Image:	



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab00119(ACTL7A antibody) at dilution of 1:800