

## ACTRT2 antibody

### Product Information


Catalog No.:	FNab00128
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	actin-related protein T2
Synonyms:	Actin-related protein T2 (ARP-T2) Actin-related protein M2 ACTRT2 ARPM2
Observed MW:	42 kDa
Uniprot ID :	Q8TDY3

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:200-1:2000
Image:	



rat testis tissue were subjected to SDS PAGE followed by western blot with FNab00128(ACTRT2 Antibody) at dilution of 1:600