

## ACY3 antibody

### Product Information

Catalog No.:	FNab00132
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Plays an important role in deacetylating mercapturic acids in kidney proximal tubules. Also acts on N-acetyl-aromatic amino acids(By similarity).

### Immunogen information

Immunogen:	aspartoacylase(aminocyclase) 3
Synonyms:	N-acyl-aromatic-L-amino acid amidohydrolase (carboxylate-forming) Acylase III Aminoacylase-3 (ACY-3) Aspartoacylase-2 Hepatitis C virus core-binding protein 1 (HCBP1, HCV core-binding protein 1) ACY3 ASPA2
Observed MW:	35 kDa
Uniprot ID :	Q96HD9

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:5000
Image:	



mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab00132(ACY3 antibody) at dilution of 1:400