

ADAMDEC1 antibody

Product Information

Catalog No.:	FNab00145
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May play an important role in the control of the immune response and during pregnancy.

Immunogen information


Immunogen:	ADAM-like, decysin 1
Synonyms:	ADAM DEC1 A disintegrin and metalloproteinase domain-like protein decysin-1 (ADAM-like protein decysin-1) ADAMDEC1
Observed MW:	53 kDa, 44 kDa
Uniprot ID :	O15204

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:1000; IHC: 1:20-1:200; IP: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using FNab00145(ADAMDEC1 Antibody) at dilution of 1:50



MCF-7 cells were subjected to SDS PAGE followed by western blot with FNab00145(ADAMDEC1 Antibody) at dilution of 1:300