

## ADAMTS12 antibody

### Product Information

Catalog No.:	FNab00146
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Metalloprotease that may play a role in the degradation of COMP. Cleaves also alpha-2 macroglobulin and aggrecan. Has anti-tumorigenic properties.

### Immunogen information


Immunogen:	ADAM metalloproteinase with thrombospondin type 1 motif, 12
Synonyms:	A disintegrin and metalloproteinase with thrombospondin motifs 12 (ADAM-TS 12, ADAM-TS12, ADAMTS-12) ADAMTS12
Observed MW:	170 kDa
Uniprot ID :	P58397

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:200-1:1000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human stomach cancer tissue slide using FNab00146(ADAMTS12 Antibody) at dilution of 1:50



mouse placenta tissue were subjected to SDS PAGE followed by western blot with FNab00146(ADAMTS12 Antibody) at dilution of 1:300