

ADCK5 antibody

Product Information

Catalog No.:	FNab00154
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


The function of this protein is not yet clear. It is not known if it has protein kinase activity and what type of substrate it would phosphorylate(Ser, Thr or Tyr).

Immunogen information

Immunogen:	aarF domain containing kinase 5
Synonyms:	Uncharacterized aarF domain-containing protein kinase 5 ADCK5
Observed MW:	66 kDa
Uniprot ID :	Q3MIX3

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:200-1:2000
Image:	



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab00154(ADCK5 antibody) at dilution of 1:800