

ADH1B antibody

Product Information

Catalog No.:	FNab00162
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	alcohol dehydrogenase 1B(class I), beta polypeptide
Synonyms:	All-trans-retinol dehydrogenase [NAD(+)] ADH1B Alcohol dehydrogenase 1B Alcohol dehydrogenase subunit beta ADH1B ADH2
Observed MW:	40 kDa
Uniprot ID :	P00325

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:5000
Image:	



mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab00162(ADH1B Antibody) at dilution of 1:1000