

ADNP antibody

Product Information

Catalog No.:	FNab00174
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Potential transcription factor. May mediate some of the neuroprotective peptide VIP-associated effects involving normal growth and cancer proliferation.

Immunogen information

Immunogen:	activity-dependent neuroprotector homeobox
Synonyms:	Activity-dependent neuroprotector homeobox protein Activity-dependent neuroprotective protein ADNP ADNP1 KIAA0784
Observed MW:	135 kDa
Uniprot ID :	Q9H2P0

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:5000; IF: 1:10-1:100
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab00174(ADNP antibody) at dilution of 1:1000