

ADRB2 antibody

Product Information

Catalog No.:	FNab00180
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Beta-adrenergic receptors mediate the catecholamine-induced activation of adenylate cyclase through the action of G proteins. The beta-2-adrenergic receptor binds epinephrine with an approximately 30-fold greater affinity than it does norepinephrine.

Immunogen information


Immunogen:	adrenergic, beta-2-, receptor, surface
Synonyms:	Beta-2 adrenergic receptor Beta-2 adrenoreceptor (Beta-2 adrenoceptor) ADRB2 ADRB2R B2AR
Observed MW:	47 kDa
Uniprot ID :	P07550

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:5000; IHC: 1:50-1:200
Image:	



Immunohistochemistry of paraffin-embedded human heart tissue slide using FNab00180(ADRB2 Antibody) at dilution of 1:50



mouse heart tissue were subjected to SDS PAGE followed by western blot with FNab00180(ADRB2 antibody) at dilution of 1:1000