

## AEBP2 antibody

### Product Information

Catalog No.:	FNab00188
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	9B3
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background

DNA-binding transcriptional repressor. May interact with and stimulate the activity of the PRC2 complex, which methylates 'Lys-9' and 'Lys-27' residues of histone H3.

### Immunogen information

Immunogen:	AE binding protein 2
Synonyms:	Zinc finger protein AEBP2 Adipocyte enhancer-binding protein 2 (AE-binding protein 2) AEBP2
Observed MW:	55 kDa
Uniprot ID :	Q6ZN18

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	

HeLa cells were subjected to SDS PAGE followed by western blot with FNab00188(AEBP2 antibody) at dilution of 1:1000

