

MLLT3 antibody

Product Information


Catalog No.:	FNab00191
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	myeloid/lymphoid or mixed-lineage leukemia(trithorax homolog, Drosophila); translocated to, 3
Synonyms:	Protein AF-9 ALL1-fused gene from chromosome 9 protein Myeloid/lymphoid or mixed-lineage leukemia translocated to chromosome 3 protein YEATS domain-containing protein 3 MLLT3 AF9 YEATS3
Observed MW:	64 kDa
Uniprot ID :	P42568

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:50-1:1000
Image:	



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab00191(MLLT3 antibody) at dilution of 1:600