

AFM antibody

Product Information

Catalog No.:	FNab00200
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Vitamin E binding protein. May transport vitamin E in body fluids under conditions where the lipoprotein system is not sufficient. May be involved in the regulation and transport of vitamin E at the blood-brain barrier.

Immunogen information


Immunogen:	afamin
Synonyms:	Afamin Alpha-albumin (Alpha-Alb) AFM ALB2 ALBA
Observed MW:	70 kDa
Uniprot ID :	P43652

Application

Reactivity:	Human
Tested Application:	ELISA, WB,IP
Recommended dilution:	WB: 1:1000-1:10000; IP: 1:1000-1:10000
Image:	



IP Result of anti-AFM (IP:FNab00200, 4ug;
Detection:FNab00200 1:2000) with human plasma
lysate 4000ug.



human blood tissue were subjected to SDS PAGE
followed by western blot with FNab00200(AFM
Antibody) at dilution of 1:2000