

AGFG2 antibody

Product Information


Catalog No.:	FNab00211
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	ArfGAP with FG repeats 2
Synonyms:	Arf-GAP domain and FG repeat-containing protein 2 HIV-1 Rev-binding protein-like protein Rev/Rex activation domain-binding protein related (RAB-R) AGFG2 HRBL RABR
Observed MW:	49 kDa
Uniprot ID :	O95081

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:200-1:2000; IF: 1:10-1:100; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lymphoma using FNab00211(AGFG2 antibody) at dilution of 1:50



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab00211(AGFG2 antibody) at dilution of 1:800