

PHYKPL antibody

Product Information

Catalog No.:	FNab00225
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	alanine-glyoxylate aminotransferase 2-like 2
Synonyms:	5-phosphohydroxy-L-lysine phospho-lyase Alanine--glyoxylate aminotransferase 2-like 2 PHYKPL AGXT2L2
Observed MW:	50 kDa, 19 kDa
Uniprot ID :	Q8IUZ5

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human pancreas using FNab00225(AGXT2L2 antibody) at dilution of 1:50