

## AHNAK antibody

### Product Information

Catalog No.:	FNab00230
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background

May be required for neuronal cell differentiation.

### Immunogen information

Immunogen:	AHNAK nucleoprotein
Synonyms:	Neuroblast differentiation-associated protein AHNAK Desmoyokin AHNAK PM227
Observed MW:	629 kDa
Uniprot ID :	Q09666

### Application

Reactivity:	Human
Tested Application:	ELISA, IF
Recommended dilution:IF:	1:50-1:500