

## AIG1 antibody

### Product Information

Catalog No.:	FNab00237
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


May play a role in androgen-regulated growth of hair follicles.

### Immunogen information


Immunogen:	androgen-induced 1
Synonyms:	Androgen-induced gene 1 protein (AIG-1) Fatty acid esters of hydroxy fatty acids hydrolase AIG1 (FAHFA hydrolase AIG1) AIG1
Observed MW:	30 kDa
Uniprot ID :	Q9NVV5

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:1000; IHC: 1:20-200
Image:	



Immunohistochemistry of paraffin-embedded human ovary tumor using FNab00237(AIG1 antibody) at dilution of 1:100



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab00237(AIG1 antibody) at dilution of 1:400